After Service Card 售後服務保障卡

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Wish you health and peace.
願健康與平安歸給您。

Dr.Rock Energy Pain Relief Patch

Product Description
Dr.Rock Energy Pain Relief Patch is a non-invasive therapeutic device that can be magnetically attached and detached from the main unit. It utilizes graphene heating technology to release far-infrared rays and negative ions. The patch also contains PEMF (Pulsed Electromagnetic Field) energy at 7.83Hz, along with a natural trace element formula in the energy patch, to achieve pain relief effects.

Six Key Features

- Graphene heating
- PEMF 7.83Hz
- 30-45°C thermostatic
- Vibration adjustment
- Magnetism
- Easy to carry

Product Diagram

Product User Manual

Loc & Func

Loc & Func	Instructions for use
Power BTN / Vibration BTN / PEMF BTN	Long press for 2 seconds to power on/off; after powering on, short press to adjust vibration intensity. Adjustment sequence: (1st gear ⏪) (2nd gear ⏩) (3rd gear ⏸) → off. After powering on, PEMF (Pulsed Electromagnetic Field) is turned on by default, with a frequency of 7.83Hz. Quick double-click to enable/disable PEMF function.
Temp Knob	Use the temperature control knob to adjust the temperature within the range of 30 to 45°C. Clockwise rotation increases the temperature by 1 degree until it reaches 45°C, then stops. Counterclockwise rotation decreases the temperature by 1 degree until it reaches 30°C, then stops. When selecting the minimum temperature, counterclockwise rotation turns off the heat therapy. When the heat therapy is turned off, rotating clockwise turns it on. "30-45°C" indicates the current heat therapy temperature setting. "00" means heat therapy is turned off. "—" means the controller is separated from the pain relief patch.
Temp Indicator light	Standard USB-C charging socket (recommended adapter: DV2A)
USB-C Port	Displaying 7.83Hz indicates that PEMF is turned on, while no display indicates that PEMF is turned off.
PEMF Indicator light	Displaying current vibration levels: (1 gear ⏪) (2 gear ⏩) (3 gear ⏸)
Vibration mode indicator light	The flashing of the "i" icon indicates that heating is in progress. When the set temperature is reached, the "i" icon disappears, indicating that the heating operation has stopped.
heating Indicator light	Display the current battery power status. Full battery: 100% to 80% - Battery indicator light is green. Normal battery level: 30% to 5% - Battery indicator light is blue. Low battery level: 30% to 5% - Battery indicator light is red. When the battery is low, there will be a "beep" sound notification. Charging status: Red breathing effect on the battery indicator light. Fully charged: Battery indicator light remains solid green.
Battery Indicator light	This product has a memory function, which retains the previous temperature setting even after power off or separation from the controller.

Product Usage Diagram

Specifications

Model: EP100
Name: Dr.Rock Energy Pain Relief Patch
Auxiliary heating: graphene membrane
Hot compress temperature: 30-45°C
Adjustable temperature: 30-45°C (Note: The above will vary by ±5°C due to environment and usage)
Vibration: 3-speed adjustment
Automatic timing: every 20 minutes. (To use again, please remove the controller and restart.)
Ultra long wave frequency: 7.83 Hz
Power input: USB-C 5V/1.2A
Lithium battery: 3.7V/1500mAh, 5.5Wh
Battery life: about 1 hour (depending on usage status)
Charging time: about 2 hours
Size: 149*79*30mm / Weight: 88g / Colour: Black

Applicable Scope

According to different painful areas, it can be applied to the knee, ankle, hand joints (De Quervain's Tenosynitis or "Mommy Thumb", Tennis Elbow, Rheumatoid Arthritis), joint pain, cervical spine, shoulder, lumbar spine, muscle soreness in the legs (post-exercise muscle soreness, musclestrain pain), and as an adjunctive treatment for hip pain.

Safety Precautions For Use

- Individuals with abnormal skin conditions, wounds, eczema, dermatitis, allergies, and pregnant women should use this product with caution.
- Discontinue use and consult a doctor if any adverse reactions occur.
- For external use only. The energy patch can be reused up to three times.
- Do not use the product if it is contaminated, if there are spherical or chunky particles falling off, or if there are any sharp or prickly sensations.

Package Contents

Dr.Rock Energy Pain Relief Patch x1	Product storage bag x1
Waist and leg pain energy patch x10	User's manual x1
Shoulder neck pain energy patch x10	Charging cable x1
Joint pain energy patch x10	

FAQs & Troubleshooting

Fault description	Reason for the anomaly	Handling measures
The product cannot power on	Check if the battery is fully charged.	Please charge the device. If it has not been used for a long time, please charge it for at least 2 hours.
The device shuts down during use	Automatic timing has ended or insufficient battery power.	Please reboot or charge again.
The product cannot charge or shows no charging indication	1. Check if the charging port is properly connected. 2. Check if the charging power source has sufficient power.	1. Please reinstall or replace the appropriate adapter. 2. Check if the charging power source is working properly.
The product doesn't vibrate	The vibration motor is not receiving power or is damaged.	For specific issues, please contact customer support for further assistance.
The product doesn't feel warm	Wait for 1-2 minutes. The heating element may not be receiving power or it could be burned out.	For specific issues, please contact customer support for further assistance.

Note: If there are any other abnormal conditions, please contact customer service promptly.

After-Sales Service

One-year warranty

From the date of purchase, the whole machine is covered by a one-year warranty, while the accompanying accessories are covered by a six-month warranty. In the event of quality issues caused by manufacturing defects that prevent normal usage, our company will be responsible for free repairs, replacements, or necessary component replacements.

Note: The following disclaimer is not applicable within this warranty scope:

- Damages caused by accidents, misuse, mishandling, or transportation;
- Unauthorized repairs of the product;
- Failure to follow the maintenance or usage instructions stated in the user manual;
- Loss exceeding the cost of the product;
- Quality issues arising from abnormal storage or protection conditions at the customer's workplace;
- Inability to provide a valid purchase receipt;
- Products purchased from unofficial or unauthorized channels;

Environmental Protection And Recycling

The product contains a built-in rechargeable battery that may contain environmentally hazardous substances. Before disposing of the product, please make sure to turn off the power and remove the battery, or take the product to an official designated product recycling center. You can also send the discarded product to our designated recycling department (please contact customer service for recycling information), where the battery will be disassembled and recycled by the repair department.

- Please do not short-circuit the battery to prevent battery combustion and fire hazards.
- Do not dispose of the battery in fire to prevent explosions or the release of toxic substances.

Safety Precautions

- Before use, carefully read the instruction manual and follow the operating instructions.
- Do not use the product alone by individuals who are unable to take care of themselves, infants, or those who are not sensitive to heat.
- When using, pay attention to the temperature after heating and select the appropriate temperature setting that is tolerable for your skin. Avoid overheating to prevent burns.
- It is strictly prohibited to use chemical solvents to clean the device body, charging cable, etc.
- Do not dismantle the device by yourself in any way; otherwise, any resulting damage will be solely your responsibility.
- If you feel the temperature is too high or if the applied area becomes red, it is recommended to lower the heat setting or temporarily discontinue use.

Maintenance & Storage

- After using the product, remove the used energy patch promptly.
- When storing the product, keep it in a dry, well-ventilated place, protected from moisture, fire, sunlight, pressure, impact, and contamination.
- Once the energy patch is removed from its packaging, it should be used immediately and should not be placed back for storage to avoid contamination and loss of adhesive properties.
- After using the main unit of the product, turn off the power and place it in a safe location, away from children's reach.
- When storing the product or not using it for an extended period, keep it in a storage bag. Avoid exposing the product to direct sunlight, as it may cause fading, deformation, or even bursting. Additionally, if not used for a long time, it is recommended to charge the main unit of the product once a month to prevent battery damage or moisture issues.
- Do not use this product for any other purposes.

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