BianStone & Graphene Far Infrared Heating Therapy SleepPad

Product Introduction

This SleepPad is made with selected natural red Bian-Stones, combined with Terahertz energy stones and a Graphene heating film for heat generation. It operates at a low voltage of 36V, ensuring safety against electric shock. Through its multi-energy formula, it improves sleep quality, promotes body recovery and regeneration, and wakes you up feeling energized every day.

Product material and function

Natural Red BianStone: Bianstone surgery, one of the six medical techniques of Chinese medicine, mentioned in many places in the "Huangdi Neijing", it is called a stone that can cure diseases, using Bianstone as a treatment tool. Representative: Chinese medicine Bian Que. After testing by modern instruments, it contains more than 40 kinds of trace elements and minerals that are beneficial to the human body, as well as more than 20 kinds of anti-aging elements such as calcium, magnesium, zinc, chromium, strontium, and selenium that are necessary for the human body.

Graphene: Graphene far-infrared rays have the functions of improving blood circulation, improving joint pain, boosting spirit, skin care and beauty.

Terahertz: It can produce a wave heat effect, increase cell vitality, adjust nerves and the body, promote body metabolism, improve the human body's magnetic field, open up the meridians in the human body, enhance the human body's mental vitality, and promote the function of sweat glands, thereby accelerating the excretion of waste in the body, beneficial to health.

Product Structure

- 1: Nylon edging
- 2: Heat-gathering layer
- 3: High-density mesh cloth
- 4: red BianStone / terahertz
- 5: Graphene heating film
- 6: Bamboo charcoal fiber cotton
- 7: 3D shockproof net

Eight characteristics:

- Natural premium BianStone
- Terahertz rough stone
- constant temperature 20~50°C
- Graphene heating film
- Far-infrared 5-15 μm
- Dual-layer full touch design
- Anti-leaking Anti-leakage
- Foldable and rollable

Product usage diagram

Product Connection Diagram:

1. Connect the six-pin connector (male) to the temperature controller's six-pin connector (female) by aligning the guide slot and following the arrow indication, then rotate it tightly.
2. Connect the other end of the temperature controller's DC connector (female) to the power supply's DC connector (male) by aligning the plug guide slot and rotating it to secure the connection.
3. Insert the matching power cord into the power interface.

Controller use and troubleshooting

Power Indicator Light, Temperature Display, Time Unit - Hours, Heating Indicator Light, Temperature Up Button, Timer Increase Button, Temperature Down Button, Timer Decrease Button, Power On/Off Button

1. After the power is turned on, long press the "ON/OFF" key to turn on or off. The default temperature "30°C" is set to "8" hours when it is used for the first time.
2. The "ON/OFF" icon is for temperature adjustment. Press "UP" or "DOWN" to adjust for temperature rise and fall. Note: After the temperature is set, the displayed temperature will automatically jump to the actual sensing temperature.
3. "ON/OFF" is for timing adjustment, press "ON/OFF" key for timing adjustment, the unit is H/Hours. Note: After the timer is counted down, it will automatically stop working and shut down when the set time is reached.
4. The icon flashes during adjustment, and the icon is always on for 6 seconds after setting to indicate completion. Note: After the temperature is set, the "ON/OFF" icon disappears to indicate that the heating is stopped.
5. The thermostat has a memory function. After the last function adjustment, the default last adjustment parameters will be used again.
6. Temperature safe mode:

mode	set up	Automatic heating time
A.Warm therapy 1	20-45°C	Can be set manually for 1-12 hours
B.Warm therapy 2	46-50°C	Heating for 4 hours, automatically return to 39°C

7. Temperature Controller Display Troubleshooting

Product Connection Diagram:

1. Connect the six-pin connector (male) to the temperature controller's six-pin connector (female) by aligning the guide slot and following the arrow indication, then rotate it tightly.
2. Connect the other end of the temperature controller's DC connector (female) to the power supply's DC connector (male) by aligning the plug guide slot and rotating it to secure the connection.
3. Insert the matching power cord into the power interface.

Controller use and troubleshooting

Power Indicator Light, Temperature Display, Time Unit - Hours, Heating Indicator Light, Temperature Up Button, Timer Increase Button, Temperature Down Button, Timer Decrease Button, Power On/Off Button

The product particularly soothes and improves the following peoples

1. Insomnia, long-term fatigue, deficiency and cold constitution;
2. Women with physical pains, cold hands and feet, edema and obesity, sub-health;
3. Back pain, shoulder and neck pain, numbness of hands and feet;
4. Low immune function, Muscle soreness, sports injuries.

Physiotherapy mode

- **Sleep Heat Therapy:** Set the temperature controller to 30-40°C. Recommended heating usage: 2-4 hours per session. Suitable for individuals with cold constitution, poor sleep quality, poor body microcirculation, chronic fatigue, rheumatic pain, cold hands and feet, and weakened immune system.
- **Warm Conditioning Therapy:** Set the temperature controller to 41-49°C. Recommended heating usage: 2-4 hours per session. Suitable for combating stubborn diseases, rheumatic pain, dysmenorrhea, cold uterus, blood stasis, cold hands and feet, poor body microcirculation, edema and obesity, and sub-healthy individuals. It helps in the recovery of symptoms such as hemorrhoids and prostatitis in men. It also aids in the recovery from sports injuries, lower back pain, bone pain, shoulder and neck pain, and numbness in the hands and feet. It can accelerate blood circulation, promote sweating and detoxification, effectively alleviate and improve pain symptoms.
- **Natural Room Therapy:** Turn off the power of the temperature controller or natural care of indoor temperature. Recommended usage: Long-term. In summer, it provides a cool and comfortable sensation. Prolonged use allows the body to absorb beneficial far-infrared and terahertz energy emitted by the BianStone, helping to eliminate cold and dampness in the body, promote metabolism, and maintain and improve the normal of the body's microcirculation system.

Caution: High blood pressure individuals and pregnant women should avoid use with caution when the temperature is between 41-50°C.

Specifications:

Model: SP719
Product name: Dr.Rock BianStone & Graphene Far Infrared Heating Therapy SleepPad
Size : L70 x W160 x H15 CM / Product net weight : 6.9kg
BianStones quantity : 205 Pcs (120 large Pcs + 85 small Pcs)
Terahertz quantity : 9 Pcs
Therapy Modes : Natural temperature therapy, Sleep therapy, Thermal therapy
Adapter input / output : 100-240V~50/60Hz, DC 36V/5.8A 208W
SleepPad input power : Maximum 208W
Current Voltage : Human body safety low voltage max 36V/5.8A
Heating device : Graphene heating film
Adjustable range of temperature controller : 20~50°C
Temperature controller can set timing : 1~12 hours
Heating temperature : 20~50°C, adjustable by lifting button
Accessories : 3D shock net, Density mesh cloth, Bamboo charcoal fiber cotton

Package Contents

- * SleepPad x1
- * Temperature controller x1
- * Storage bag x1
- * Manual book x1
- * Charger adapter x1
- * Q&A Card x1

Warranty description

Dr.Rock provides a quality warranty service for hardware products, accessory materials, and craftsmanship defects in the packaging. The warranty period for the SleepPad Pro (including the temperature controller) is one year from the date of purchase. Normal wear and tear, aging, as well as damages caused by accidents, improper modification, or misuse, are not covered under this warranty. To obtain after-sales service, please contact the authorized brand distributor of Dr.Rock, visit an authorized retailer, or contact customer service online, while providing relevant purchase proof. (invoice or purchase record)

Product maintenance methods

You can use a wet towel to wipe the contaminated part when it is not powered on, or use a disinfectant wet towel to disinfect it regularly, and then turn it on for heating after drying.

Precautions

1. Read the manual carefully before use and use it in accordance with the operating instructions.
2. Don't let it be used alone by people who cannot take care of themselves, infants or children who are not sensitive to heat.
3. It is forbidden to use in water and humid environment.
4. Do not dry clean or wash. You can wipe the contaminated part with a wet towel when the power is not turned on, and then turn on the power to heat it after drying.
5. Do not insert sharp objects into the product.
6. Do not place bedding and other items in a small area to cover the heating area when the product is heated to prevent local overheating and causing malfunctions.
7. It is forbidden to place flammable materials such as fireworks and lighters on the surface of the product to prevent fire.
8. Please confirm whether the power supply is cut off before going out.
9. Using the heating function, the temperature can be adjusted according to the tolerance of the skin. In order to prevent skin burns, if you feel unwell, stop using it in time.

Disclaimer

1. The use of a heating pad carries inherent risks, and users should follow the operating instructions and maintenance guidelines provided in the product manual. We are not responsible for any issues resulting from user error or failure to follow the instructions, and users assume all risks associated with the use of the heating pad.
2. The heating pad may not be suitable for certain individuals, including infants, individuals with head tumors, leukemia patients, individuals with bleeding disorders, individuals with metal implants or pacemakers, individuals with severe cardiopulmonary insufficiency, individuals with epilepsy or mental disorders, individuals with undiagnosed acute abdominal pain, pregnant women, women in the gestational or lactation period, individuals with heat insensitivity, or individuals with other specific health conditions. It is advised to consult a doctor's advice or use the heating pad under the guidance of a doctor.
3. We make no express or implied warranties regarding the performance, suitability, quality, or reliability of the heating pad. We are not liable for any direct or indirect loss or damage resulting from the use or inability to use the heating pad, nor are we liable for any personal injury, property damage, or other adverse consequences arising from the use of the heating pad.
4. When using the heating pad, exercise caution and remain vigilant for any abnormal conditions such as excessive heat, overheating of the power cord, or other safety issues. If any issues are detected, discontinue use immediately and seek professional assistance.

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