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
### BianStone & Graphene Far Infrared Heating Therapy SleepPad Pro Ten characteristics Product Introduction THz This SleepPad Pro is made with selected natural red Bian-Stones, combined with erahertz energy stones and a Graphene heating film for heat generation, as well as PEMF (Pulsed Electromagnetic Field) 7.83Hz Schumann wave energy. It operates at a low PEME voltage of 36V, ensuring safety against electric shock. Through its multi-energy formula, 7.83 Hz it improves sleep quality, promotes body recovery and regeneration, and wakes you up Product material and function Dual-layer full touch design Natural Red BianStone: Bianstone surgery, one of the six medical techniques of Chinese medicine, mentioned in many places in the "Huangdi Neijing", it is called a stone that can cure diseases, using Bianstone as a treatment tool Representative: Chinese medicine Bian Que. After testing by modern instruments, it contains more than 40 kinds of trace elements and minerals that are beneficial to the human body, as well as more than 20 Product usage diagram kinds of anti-aging elements such as calcium, magnesium, zinc, chromium, strontium, and selenium that are necessary for the human body. Graphene: Nobel Prize in Physics, a revolutionary material for biomedicine in the 21st century. Graphene far-infrared rays have the functions of improving blood circulation, **PEMF**: World Class Energy Technology. The cells of the human body need energy to achieve optimal functional status. Pulsed magnetic fields are a soothing and energetic method that can help the human body replenish energy, activate the vitality of the body, and dredge the meridians. Terahertz: It can produce a wave heat effect, increase cell vitality, adjust nerves and the body, promote body metabolism, improve the human body's magnetic field, open up the meridians in the human body, enhance the human body's mental vitality, and promote the function of sweat glands, thereby accelerating the excretion of waste in the body, beneficial to health. connection port 3 Product Structure 1: Nylon edging 3: High-density mesh cloth → 5: Natural red BianStone / terahertz 6: Graphene heating film 7: Bamboo charcoal fiber cotton layer

### **Product Connection Diagram:**

- 1 Connect the six-pin connector (male) to the temperature controller's six-pin connector (female) by aligning the guide slot and following the arrow indication, then rotate it tightly.
- (2) Connect the other end of the temperature controller's DC connector (female) to the power supply's DC connector (male) by aligning the plug guide slot and rotating it to secure the connection. 3 Insert the matching power cord into the power interface.

### Thermostat Usage and Troubleshooting

- 2Hz Deep Sleep Mode Time Unit - Hours Heating Indicator Light

1. After the power is turned on, long press the "(b)" key to turn on or off. The default temperature "30°C" is set to "8" hours when it is used for the first time; 2. The " icon is for temperature adjustment. Press " or " or " or " to adjust for temperature rise and fall. Note: After the temperature is set, the displayed temperature will automatically jump to the actual sensing temperature; 3. "  $\overset{ullet}{ullet}$ " is for timing adjustment, press " $ilde{ riangle}$ " or " $ilde{ riangle}$ " key for timing adjustment

the unit is H/hour. Note: After the timer is counted down, it will automatically stop working and shut down when the set time is reached; 4. The icon flashes during adjustment, and the icon is always on for 6 seconds after setting to indicate completion. Note: After the temperature is set, the " — " icon flashes to indicate heating, and when the set temperature is reached, the " ... "

icon disappears to indicate that the heating is stopped;
5. PEMF regulation: default off. Press PEMF button to turn on. After 15 minutes in 7.83Hz sleep mode, it will automatically switch to the 2Hz deep sleep mode and continue running for another 15 minutes before automatically ending.

6. The thermostat has a memory function. After the last function adjustment, the default last adjustment parameters will be used again; 7. Temperature mode set up Automatic heating time safe mode: A.Warm therapy 1 20-45°C Can be set manually for 1-12 hours

B.Warm therapy 2 46–50°C Heating for 4 hours, automatically return to 39°C

8. Temperature Controller Display Troubleshooting

Display	Cause of issue		Troubleshooting
E1	The thermostat is not connected to the product or the connection is in poorcontact;     The temperature sensor probe fails.		Check the connection line, whether the connection is normal, and it needs to be re-energized.     If the product is still abnormal, please contact after-sales service.
EH	High temperature error or temperature out of protection		Turn off the power, turn on the product again after cooling to confirm, if it shows that the product is abnormal, please contact after-sales
EU	Abnormal voltage	Check the power supply, whether it is normal;     Recommended output power ≥ 208W certified power adapter or car charger;     If there is any abnormality after reinserting the power supply and restarting, please contact the after–sales service.	
EC	Current overcurrent		

# Physiotherapy mode

- Sleep Heat Therapy: Set the temperature controller to 30-40°C. If necessary, turn on the PEMF switch (default 7.83Hz sleep mode, automatically switches to PHz deep sleep mode after 15 minutes, continues for another 15 minutes, and automatically ends). Recommended heating usage: 6-12 hours per session. Suitable for individuals with cold constitution, poor sleep quality, poor body microcirculation, chronic fatigue, rheumatic pain, cold hands and feet, and weakened immune system.
- Warm Conditioning Therapy: Set the temperature controller to 41–49°C. If necessary, turn on the PEMF switch (default 7.83Hz sleep mode, automatically switches to 2Hz deep sleep mode after 15 minutes, continues for another 15 minutes, and automatically ends). Recommended heating usage: 2-4 hours per session. Suitable for combating stubborn diseases, rheumatic pain, dysmenorrhea, cold uterus, blood stasis, cold hands and feet, poor body microcirculation, edema and obesity, and sub-healthy individuals. It helps in the recovery of symptoms such as hemorrhoids and prostatitis in men. It also aids in the recovery from sports injuries, lower back pain, bone pain, shoulder and neck pain, and numbness in the hands and feet. It can accelerate blood circulation, promote
- sweating and detoxification, effectively alleviate and improve pain symptoms. Natural Room Therapy: Turn off the power of the temperature controller (maintain room temperature during summer) or natural care of indoor temperature. Recommended usage: Long-term.
- In summer, it provides a cool and comfortable sensation. Prolonged use allows the body to absorb beneficial far-infrared and terahertz energy emitted by the BianStone, helping to eliminate cold and dampness in the body, promote metabolism and maintain and improve the normal of the body's microcirculation system.
- Caution: High blood pressure individuals and pregnant women should avoid or use with caution when the temperature is between 41-50°C.

### Specifications

Model: SP716

Product name: Dr.Rock BianStone & Graphene Far Infrared Heating Therapy SleepPad Pro Size: L70 x W160 x H1.5 CM Product net weight: 8.8kg

BianStones quantity: 253 Pcs (166 large Pcs + 87 small Pcs ) Terahertz quantity: 50 Pcs PEMF (Pulsed Magnetic Field) transmitter: 1 Group

Therapy Modes: Natural temperature therapy, Sleep therapy, Thermal therapy, PEMF therapy. Adapter input/output: 100-240V~50/60Hz, DC 36V/5.8A 208W SleepPad Pro input power: Maximum 208W Current Voltage: Human body safety low voltage max 36V/5.8A

PEMF adjustable gear: 7.83Hz, 2Hz Heating device: Graphene heating film Adjustable range of temperature controller: 20~50°C Temperature controller can set timing: 1~12 hours Heating temperature: 20~50°C, adjustable by lifting button Accessories: 3D shock net, Density mesh cloth,

# Bamboo charcoal fiber cotton (fireproof) The product particularly soothes and improves the following peoples

- Insomnia, chronic fatigue, and cold body constitution. 2. Female physiological discomfort, cold extremities, edema, obesity.
- 3. Low immune function, suboptimal health. 4. Lower back pain, shoulder and neck pain, numbness in hands and feet.
- 5. Muscle soreness, sports injuries.

# Package Contents

\* SleepPad Pro x1 \* Temperature controller x1 \* Storage bag x1 \* Charger adapter x1\* Manual book x1 \* Q&A Card x1

Warranty description

Dr.Rock provides a quality warranty service for hardware products, accessory materials, and craftsmanship defects in the packaging. The warranty period for the SleepPad Pro (including the temperature controller) is one year from the date of purchase. Normal wear and tear, aging, as well as damages caused by accidents, improper modification, or misuse, are not covered under this warranty. To obtain after-sales service, please contact the authorized brand distributor of Dr.Rock, visit an authorized retailer, or contact customer service online, while providing relevant purchase proof.

(invoice or purchase record)

## **Product Maintenance and Care Methods**

You can clean the contaminated parts with a damp cloth when the power is not connected, or regularly disinfect them with disinfectant wipes. After cleaning, let it air dry before reconnecting the power and using it for heating.

# A Precautions

- 1. Carefully read the instruction manual before use and follow the operating
- 2. Do not use the product alone by individuals who are unable to care for themselves, infants, or those who are not sensitive to heat.
- Do not use in wet or damp environments 4. Do not dry clean or wash. You can clean the contaminated parts with a damp
- cloth when the power is not connected. Let it air dry before reconnecting the power and using it for heating.
- 5. Do not insert sharp objects into the product.

stop using it immediately.

- 6. Do not place bedding or other items that cover only a small area of the heating zone while the product is heating to prevent local overheating and malfunction.
- 7. Do not place flammable items such as fireworks, firecrackers, or lighters on
- the surface of the product to prevent fire hazards. 8. Before finishing use or leaving the area, ensure that the power is disconnected. 9. When using the heating function, adjust the temperature according to your skin tolerance level. To prevent burns, adjust the temperature settings based on your individual body sensitivity. If you experience any discomfort,

### a doctor's advice or use the heating pad under the guidance of a doctor. 3. We make no express or implied warranties regarding the performance suitability, quality, or reliability of the heating pad. We are not liable for any

the use of the heating pad.

Disclaimer

the heating pad, nor are we liable for any personal injury, property damage, or other adverse consequences arising from the use of the heating pad. 4. When using the heating pad, exercise caution and remain vigilant for any

1. The use of a heating pad carries inherent risks, and users should follow the

manual. We are not responsible for any issues resulting from user error or

2. The heating pad may not be suitable for certain individuals, including infants,

individuals with head tumors, leukemia patients, individuals with bleeding

disorders, individuals with metal implants or pacemakers, individuals with

severe cardiopulmonary insufficiency, individuals with epilepsy or mental

disorders, individuals with undiagnosed acute abdominal pain, pregnant women,

women in the gestational or lactation period, individuals with heat insensitivity

or individuals with other specific health conditions. It is advised to consult

direct or indirect loss or damage resulting from the use or inability to use

failure to follow the instructions, and users assume all risks associated with

operating instructions and maintenance guidelines provided in the product

abnormal conditions such as excessive heat, overheating of the power cord, or other safety issues. If any issues are detected, discontinue use immediately and seek professional assistance.

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